




SPECIFICATION		FEATURES
POWER SUPPLY	220-240VAC 50/60Hz 20W IEC Fused mains inlet.	Integrated peristaltic pump with Tygor
UV LAMP	185nm/254nm Wavelength, 4 Watt	Countdown timer with Alarm
SAMPLE VOLUME	50 - 250 ml	IP65 Safety interlocked enclosure
SAMPLE TYPE	Potable, surface and ground water	Variable sample volume
INTENDED USE	Digestion of water samples for subsequent detection of Heavy metal content (Not suitable for water disinfection)	UV lamp indicator
		Replaceable lamp
SAMPLE PUMP	Peristaltic pump - 12ml/min flow rate	
SAMPLING TUBE	4mm O/D 2mm ID Tygon tubing	

Supplied with: 5 Sample beakers, UK IEC Mains lead, 0.5m Tygon tubing

## **APPLICATION**

Surface and ground waters can contain complex organic compounds from vegetation, animals, biological effluents, soils etc. These organic compounds form bonds with metal ions which can later be released into the food chain. The Metalyser series of instruments measure the metal content in the water, however, it is not possible to measure the organically bound metal content.

Figure 1 shows a voltammogram for a typical river sample containing low levels of Cadmium and Lead. The hump in the graph indicates the presence of the metals but is not quantifiable.






FIGURE 1

FIGURE 2

Trace20 has developed the HMUV1825 digester to break down the organic compounds and in doing so release the heavy metals back into the water. Figure 2 shows the voltammogram for the sample after 30 minutes digestion. The hump has now separated into two sharp peaks for each of the metals which are quantifiable.

## ORDERING INFORMATION

UV Digester HMUV1825 - Laboratory

Aqueous sample digester designed to break down highly organic ground and surface water samples, for subsequent analysis with laboratory Metalysers. 220-240V

➤ UV Digester HMUV1825P - Portable

Aqueous sample digester designed to break down highly organic ground and surface water samples, for subsequent analysis with laboratory Metalysers. Battery powered.



**HMUV1825** 

n tubing



HMUV1825P

Trace<sub>e</sub>o Limited Technology Centre, Unit 4, Transigo, Gables Way, Thatcham, Berkshire, RG19 4ZA, United Kingdom

> Tel +44 (0) 1635 866772 Fax +44 (0) 1635 873509 sales@trace2o.com www.trace2o.com



